

HENROGRIP®



USER MANUAL

HENROGRIP® FROM THE AGE OF 8 YEARS OLD

Oktober 2023

Table of content

Important information	3
Intended use	3
Storage	3
Instruction for use.....	4
Precautions.....	4
Daily operation	4
Warnings	4
Environmentally friendly disposal	4
Technical datasheet	5

Important information

Before using HENROGRIP® read these instructions closely.

Intended use

HENROGRIP® is designed and intended for food only.

HENROGRIP® has been developed by Danish occupational therapists to ensure a high-quality cutlery for people that experience loss of strength and motoric skills in eating situations. HENDYGRIP® Is BPA Free and can safely be used by children under the age of 3 years.

HENROGRIP® is designed with a curved grip and a material that gives a firm grip. The cutlery have been designed for both left and right-handed.

Storage

Store dry, at room temperature and out of sunlight.

Instruction for use

1. Unpack the HENROGRIP® and dispose of the packing properly.
2. Before use wash HENROGRIP®.
3. The HENROGRIP® are now ready for use.

Precautions

After unpacking the grip always wash the HENROGRIP® before use.

The grip can be used for all foods.

Do not heat the grip over a candle light or fire.

Daily operation

HENROGRIP® can be washed in Dish Washer with up to 65° C degrees.

Warnings



PLEASE DO NOTICE WASHING THE HENROGRIP IN INDUSTRIAL DISHWASHER THAT WASHES ON 90° C, MAY DAMAGES THE HANDLE ON THE GRIP.

Environmentally friendly disposal

Please dispose the cup in an environmentally friendly neutral way.



HENROGRIP® are made of material which should be recycled. Please therefore dispose the old HENROGRIP® via suitable collection systems.



IF YOU HAVE ANY DOUBT OF WHEN TO DISPOSE THE HENROGRIP®, YOU CAN SEND AN EMAIL TO HENROTEK@HENROTEK.COM

Technical datasheet	
Material	PP (Polypropylene with 30 % minerals) TPE (Thermoplastic Elastomer)
BPA	HENROGRIP® is BPA Free
Marked	CE MD
Weight	40,5 G (fork) 41 G (knife) 42,5 G (Spoon)
Length	16,5 cm
Color	Grey, Red
Dishwasher	65° C degrees
Basic UDI-DI No.	57048313BY
UDI-DI HenroGrip Red Knife	05704831128319
UDI-DI HenroGrip Grey Knife	05704831128326
UDI-DI HenroGrip Grey Fork Left	05704831128340
UDI-DI HenroGrip Grey Spoon Right	05704831128364
UDI-DI HenroGrip Grey Spoon Left	05704831128388
UDI-DI HenroGrip Grey Fork Right	05704831128401
UDI-DI HenroGrip Grey Fork Left	05704831128432
UDI-DI Red Spoon Right	05704831128449
UDI-DI Red Spoon Left	05704831128456
Manufacture	Henrotek A/S, Vestermarksvej 5, 6630 Roedding

We reserve the right to make changes to specifications without notice. Subject to confirmation, availability and errors.

Standards and Regulations	<p>Regulation (EC) 1935/2004 of 27. October 2004</p> <p>Regulation (EC) 2023/2006 of 22 December 2006</p> <p>Danish executive order, BEK nr 681 of 25/05/2020</p> <p>Regulation (EC) 10/2011 of 14 January 2011</p> <p>Regulation (EU) 2017/745 of April 2017</p>
----------------------------------	---